

Global settings (administration)

This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

Global settings has 6 fields:

1. Language
2. Application IconNo file chosen
3. Email For Help
4. Conference Server Url
5. Snapshot Filename Template

Annotation settings has 5 fields:

1. Show Annotation Dialog
2. Auto Save Annotations
3. Restart Annotations Drawing
4. Show Slide Info Button
5. Show Full Slide Path In Viewport

Edit Global Settings

Global application settings

Language: en (English)

Show Research Use Only Dialog:

Application Icon: Choose File No file chosen

Email For Help: help@pathomation.com

Conference Server Url: https://test.pathomation.com/pma/test

Snapshot Filename Template: snapshot_{date}_{time}

Annotation-related settings

Show Annotation Dialog:

Auto Save Annotations:

Restart Annotations Drawing:

Show Slide Info Button:

Show Full Slide Path In Viewport:

Save Cancel

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_global&rev=1680811111

Last update: **2023/04/06 22:58**

