

Shared links (administration)

This page list down all the shared links created for a slide or slide folder. It provides detail about each of the shared links that include link type, is password protected and expiry date. You can also delete or suspend a shared link.

PATHOMATION								Sign
Shared Links								
Id	Name	Comment	Type	Is protected by password	Expires At	Actions		
03UR6dme			RCI	Not Protected	Jan 01 0001 00:00 Expired			
04qVAKU8	testje Cleg		Slide	Not Protected	Never			
05mKU1UQ			Slide	Not Protected	Never			
0UGVFCU			Grid	Not Protected	Never			
0JlMrT			Slide	Not Protected	Never			
0SHCT7aY	Cleg test		RCI	Protected	Jul 08 2021 21:00 Expired			
0wZ6CH9			RCI	Not Protected	Never			
1AqpsSq			RCI	Not Protected	Never			
113GdHMg			Slide	Not Protected	Aug 17 2020 19:00 Expired			
1Nc933W			Slide	Not Protected	Never			
1KPhw81E			Grid	Not Protected	Never			
1M3F0My			Slide	Not Protected	Never			

From: <https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link: https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_links

Last update: **2023/04/06 22:57**

