


# Shared links (administration)

This page list down all the shared links created for a slide or slide folder. It provides detail about each of the shared links that include link type, is password protected and expiry date. You can also delete or suspend a shared link.



| ID         | Name       | Summary | Type | Is password protected | Expires at        | Actions            |
|------------|------------|---------|------|-----------------------|-------------------|--------------------|
| 1000000001 | 1000000001 |         | URL  | Not Protected         | Jul 31 2025 00:00 | [Suspend] [Delete] |
| 1000000002 | 1000000002 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000003 | 1000000003 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000004 | 1000000004 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000005 | 1000000005 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000006 | 1000000006 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000007 | 1000000007 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000008 | 1000000008 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000009 | 1000000009 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000010 | 1000000010 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000011 | 1000000011 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000012 | 1000000012 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000013 | 1000000013 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000014 | 1000000014 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000015 | 1000000015 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000016 | 1000000016 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000017 | 1000000017 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000018 | 1000000018 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000019 | 1000000019 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |
| 1000000020 | 1000000020 |         | URL  | Not Protected         | None              | [Suspend] [Delete] |

From:  
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin\\_links&rev=1648661905](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_links&rev=1648661905)

Last update: **2022/03/30 20:38**

