


Ribbon layout (administration)

This admin page allows you to setup default Ribbon layout that consists of different tabs and controls. Each tab consists of different controls associated to it and you can configure appearance of tabs and controls using this page.



Below are some of the tabs and corresponding controls used in PMA.Studio.

Home tab has 7 sections:

- 1. View
- 2. Grid
- 3. Directory
- 4. Slide
- 5. Share
- 6. Search
- 7. Conferencing



Annotations tab (available in PMA.Studio) supplies instruments for working with annotations and has 3 sections:

- 1. Annotations
- 2. Shapes
- 3. Export

Grid tab supplies instruments to manage grid and it has 6 sections:

- 1. ItemSize
- 2. Synchronize
- 3. Rows/Columns
- 4. Slide
- 5. Snapshot
- 6. Share

Configure tab has 2 sections:

1. Server (available in PMA.Studio)
2. Layout

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:

https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_ribbon&rev=1648658240

Last update: **2022/03/30 19:37**

