


# Ribbon layout (administration)

This admin page allows you to setup default Ribbon layout that consists of different tabs and controls. Each tab consists of different controls associated to it and you can configure appearance of tabs and controls using this page.



Below are some of the tabs and corresponding controls used in PMA.Studio.

**Home tab** has 7 sections:

1. View
2. Grid
3. Directory
4. Slide
5. Share
6. Search
7. Conferencing



**Annotations tab (available in PMA.Studio)** supplies instruments for working with annotations and has 3 sections:

1. Annotations
2. Shapes
3. Export

**Grid tab** supplies instruments to manage grid and it has 6 sections:

1. ItemSize
2. Synchronize
3. Rows/Columns
4. Slide
5. Snapshot
6. Share

**Configure tab** has 2 sections:

1. Server (available in PMA.Studio)
2. Layout

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin\\_ribbon&rev=1648658240](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_ribbon&rev=1648658240)

Last update: **2022/03/30 19:37**

