

# Connectivity

We live in a connected world, where vertical data-silos unfortunately are still much to frequent.

With our PMA.vue and PMA.studio products, we want to help you integrate disparate information sources, and bring them together in a single harmonized interface.

## PMA.vue

With PMA.vue, you can integrate virtual slides ([find out what makes these data so unique from other imaging data](#) from [any number of slide scanners](#) into a single user interface.

You get tree-based navigation on top of your slide repositories (organized through a single PMA.core back-end [tile server](#)).

Not only can you connect to data within your own organization, PMA.vue also makes it easy to reach out to external parties and connect other to *you* in return. PMA.vue offers many sharing features, and comes with built-in snapshot capabilities to prepare manuscripts and other publications.

## PMA.studio

PMA.studio adds connectivity where PMA.vue stops.

PMA.studio allows interaction with multiple PMA.core tile servers simultaneously (as well as [PMA.start](#) and [My Pathomation](#)).

## Compare PMA.vue and PMA.studio

From:  
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=connectivity&rev=1671015132>

Last update: **2022/12/14 13:52**

