


Context menu

The product indication label in the top-left corner of the PMA.studio and PMA.vue user interface serves as a drop-down menu.



The items in the menu are the same for PMA.studio and PMA.vue.

About us

Show EULA

Manual

Keyboard shortcuts

File formats

Helpdesk

Change password

Edit settings

Sign out

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=context_menu&rev=1671024035

Last update: **2022/12/14 16:20**

