


# The Home Tab

The Home tab contains controls for:

- Layouts
- Viewport controls
- A small number of Annotations features
- Folder tools
- Grid tools
- Snapshot tools
- Magnification controls
- Tray
- Sharelinks



From:  
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=home&rev=1731502992>

Last update: **2024/11/13 16:03**

