


The Home Tab

The Home tab contains controls for:

- [Layouts](#)
- [Viewport controls](#)
- A small number of Annotations features
- [Folder tools](#)
- [Grid tools](#)
- [Snapshot tools](#)
- Magnification controls
- [Tray](#)
- [Sharelinks](#)



From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=home&rev=1733487197>

Last update: **2024/12/06 15:13**

