


HTTP vs HTTPS

PMA.studio and PMA.vue preferably are run under the [HTTPS-protocol](#) (provided you have an [SSL certificate](#) enabled on your webserver).

If you accidentally launch either application under the HTTP protocol instead of the HTTPS protocol, a notification shows during the login procedure.



The notification appears whether you intended to launch the application under HTTP or whether it was an accident. It's up to the end-user to decide whether the message is justified or not (in which case you can continue using the software as is, e.g. when running it on *localhost*).

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=https&rev=1671023463>

Last update: **2022/12/14 16:11**

