

# Keyboard Shortcuts

Below you can find a list of Keyboard Shortcuts.

Quick-Connect to server

- ctrl + q

Save snapshot

- ctrl + s

Share

- ctrl + shift + s

ctrl + m Toggle multi-selection mode [ Decrease zoom ] Increase zoom ctrl + Select next image ctrl + Select previous image ctrl + shift + Select next directory ctrl + shift + Select previous directory ctrl + shift + a Export Annotations shift + drag zoom Drag zoom alt + rotate Rotate slide

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=keyboard\\_shortcuts&rev=1734356354](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=keyboard_shortcuts&rev=1734356354)

Last update: **2024/12/16 16:39**

