


# Bookmark Panel

Bookmark panel allows user to bookmark slides. From the bookmark panel user can click on the “+” icon to add a slide to bookmark and panel lists down bookmarks.



From:  
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=panel\\_bookmark&rev=1746972921](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=panel_bookmark&rev=1746972921)

Last update: **2025/05/11 17:15**

