

Whats new?

PIMS LS 3.0.1

Added further security to the API endpoints

- We ran certain security scans and made our API endpoints lot more secure

Bug fixes

- Synchronization interoperability between “Caption” and “Classification” values in the Annotations panel and text lines pop-up window
- Disappearance of the annotations' opacity settings
- Synchronization of “Outline Color Opacity” for existing and new annotations
- Synchronization between the “Show/Hide Annotations” in the Annotation panel and the “Annotations” toggle button in the Viewport drop-down menu
- The viewport drop-down menu buttons are not working correctly
- Change slide path to slide name
- Pop out selected slide doesn't function correctly
- Add overlay grid to PIMSlS

Improvements to existing features

- Added user account information to the upper right corner along with email address and moved PIMS LS brand title to the top left
- Also moved slide loading messages and indicator to the top center of the screen, instead of top right

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=release_notes&rev=1751634696

Last update: **2025/07/04 16:11**

