

Whats new?

PIMS LS 3.0.2 Patch

Added option to allow importing PNG and TIFF images

- Added functionality to allow PNG and TIFF images annotation import. Previously application was handling well know text(.wkt), visio pharam (.mld), indica labs (.halo) and .csv formats.

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:

https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=release_notes&rev=1759844001

Last update: **2025/10/07 16:33**

