

## Whats new?

### PIMS LS 3.0.2 Patch

#### Added option to allow importing PNG and TIFF images

- Added support for importing annotations from PNG and TIFF image formats.
- Previously, the application supported annotation imports only from Well-Known Text (.wkt), Visio Pharam (.mld), Indica Labs (.halo), and .csv formats.

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=release\\_notes&rev=1759845241](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=release_notes&rev=1759845241)

Last update: **2025/10/07 16:54**

