

# Ribbon

Ribbon is a panel that sticks to the top of the screen and consists of tabs and controls. Appearance of tabs and controls are being controlled by ribbon layout, which can be specified in the Admin area.

## File Menu tab


This is the starting tab of PMA.studio and provides set of menu items used by the application. It includes:

1. About us
2. Show EULA
3. Manual
4. Keyboard shortcuts
5. File formats

## HOME tab

HOME tab consists of different sections that contain set of controls for:

1. Layout
2. Viewport
3. Navigate
4. Grid
5. Snap
6. Magnify
7. Slide
8. Tray
9. Share



## GRID tab

GRID tab consists of different sections that contain set of controls:

1. Grid Size
2. Synchronize
3. Rows and Columns (R/C)
4. Slide Zoom in and Zoom out
5. Snapshot
6. Share



## ANNOTATIONS tab

ANNOTATIONS tab consists of different sections that contain set of controls for:

1. Annotations
2. Style

- 3. Select Annotation
- 4. Simple Shapes
- 5. Polygon Shapes
- 6. Export
- 7. Import



### PRE-SET ANNOTATIONS tab

PRE-SET ANNOTATIONS tab consists of different sections that contain set of controls for:

- 1. Annotations
- 2. Stamps
- 3. Protocols



### CONFERRENCING tab

ANNOTATIONS tab consists of different sections that contain set of controls for:

- 1. Conferencing
- 2. Broadcasting



### CONFIGURE tab

CONFIGURE tab consists of different sections that contain set of controls for:

- 1. Server
- 2. Panels
- 3. Pre-set layouts
- 4. Viewport



From: <https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link: <https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon&rev=1651177715>

Last update: **2022/04/28 23:28**

