

Ribbon > Configure tab > Panels group > Toggle conference panel

Button identification

Associated command

panel-toggle-conference

Function description

Brings up the [conference panel](#):



This panel contains [experimental features](#).

Related functionality

[Ribbon > Home tab](#)

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_config_panels_toggle_conference&rev=1677849004

Last update: **2023/03/03 16:10**

