

Ribbon > Configure tab > Panels group > Toggle shared links panel

Button identification

Associated command

none

Function description

Brings up the [shared links panel](#):

Shared Links ×

Remove selected

Revoke selected

Search here...

☐	Id	Name	Comment	Type	Is protected by password	Expires At	Actions
☐	05mKU1UQ			Slide	Not Protected	Never	      
☐	2KFQwr9s			Slide	Not Protected	Never	      
☐	AcYsZzJW			Directory	Not Protected	Never	      
☐	f4PktLHb			Slide	Not Protected	Never	      

Related functionality

[Ribbon > Home tab](#)

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_config_panels_toggle_shared_links

Last update: 2023/02/27 15:10

