

Ribbon > Home tab > Layout group > Default layout

Button identification




Associated command

simple-layout

Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on a left an a [viewport panel](#) on the right.



Related functionality

[Ribbon > Home tab](#)

Last update: 2023/02/10 16:31 ribbon_home_layout_default https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_home_layout_default&rev=1676035914

From: <https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link: https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_home_layout_default&rev=1676035914

Last update: **2023/02/10 16:31**

