

# Ribbon > Home tab > Viewport group > Toggle slide label

## Button identification



## Associated command

viewport-toggle-label

## Function description

This function shows or hides the full path reference to the current slide in the top-right area of the ribbon

When enabled, the ribbon looks like this:



When disabled, the widget is hidden from the ribbon:



Unlike other toggle functionality, this one does not impact the viewport area

## Related functionality

[Ribbon > Home tab](#)

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon\\_home\\_viewport\\_toggle\\_label&rev=1676050515](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_home_viewport_toggle_label&rev=1676050515)

Last update: 2023/02/10 20:35

