

Ribbon > Home tab > Viewport group > Reset panel layout

Button identification




Associated command

reset - layout


Function description

This resets the panel layout to the default settings, as defined in the [admin interface](#).

When your layout is like this e.g.:



After clicking the button the organization of the panels resets to:



The currently loaded slides are kept.

Also note that [the configuration of the panels as stipulated through the admin interface](#) can be different from the one defined by the [Default layout button](#).

Related functionality

[Ribbon > Home tab](#)

[Default layout button](#)

[Administrative interface to define default panel configuration](#)

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_home_viewport_toggle_reset_panel_layout

Last update: 2023/02/10 18:40

