

Ribbon > Home tab > Viewport group > Toggle zoom controls

Button identification




Associated command


viewport-toggle-zoom

Function description

When enabled, the viewport looks like this:



When disabled, the widget is hidden from view:



Related functionality

[Ribbon > Home tab](#)

From:
<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=ribbon_home_viewport_toggle_zoom

Last update: 2023/02/10 19:20

