



Rotation Tools

Rotation dropdown consists of different rotation options for the viewport including flipping viewport horizontally, vertically, rotating it clockwise or counter clockwise and then reset the flips and rotations.



- Flip horizontally - flips viewport slide horizontally
- Flip vertically - flips viewport slide vertically
- Rotate counter-clockwise 90 degree
- Rotate 180 degree
- Reset flip
- Reset rotation



From:
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:
<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=rotation&rev=1747059890>

Last update: **2025/05/12 17:24**

