


# Grid: Slide Tool Group

The Slide dropdown contains tools for turning changing brightness/contrast of your current image, rotating and magnifying the slide.



Clear current cell - clears current cell of the grid  
Clear all grid cells - clears all cells of a grid  
Open slide from current cell into a new viewport panel  
Popout selected slide in a new window

From:  
<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=slide&rev=1746905182>

Last update: **2025/05/10 22:26**

