


Grid: Slide Tool Group

The Slide dropdown contains tools for turning changing brightness/contrast of your current image, rotating and magnifying the slide.



Clear current cell - clears current cell of the grid
Clear all grid cells - clears all cells of a grid
Open slide from current cell into a new viewport panel
Popout selected slide in a new window

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:

<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=slide&rev=1746905182>

Last update: **2025/05/10 22:26**

