

# Users (administration)

This page lists down all the PMA.studio users. Records can be filtered through PMA.core server url, username and date of visit. Each row consists of:

1. PMA.Studio username
2. PMA.core server url
3. PMA.core credentials used to connect with PMA.core server
4. First ever login to that PMA.core server
5. Most recent login to that PMA.core server

PMA Studio User	Server url	PMA core Credentials used	First Login	Most Recent Login
admin	https://pma.studio.pathomation.com/	admin	Monday, March 28, 2022, 10:58:56 pm	Monday, March 27, 2022, 10:51:44 pm
u_admin1	https://pma.studio.pathomation.com/	u_admin1	Monday, March 27, 2022, 10:58:56 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin2	https://pma.studio.pathomation.com/	u_admin2	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin3	https://pma.studio.pathomation.com/	u_admin3	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin4	https://pma.studio.pathomation.com/	u_admin4	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin5	https://pma.studio.pathomation.com/	u_admin5	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin6	https://pma.studio.pathomation.com/	u_admin6	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin7	https://pma.studio.pathomation.com/	u_admin7	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin8	https://pma.studio.pathomation.com/	u_admin8	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin9	https://pma.studio.pathomation.com/	u_admin9	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin10	https://pma.studio.pathomation.com/	u_admin10	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin11	https://pma.studio.pathomation.com/	u_admin11	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin12	https://pma.studio.pathomation.com/	u_admin12	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin13	https://pma.studio.pathomation.com/	u_admin13	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin14	https://pma.studio.pathomation.com/	u_admin14	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin15	https://pma.studio.pathomation.com/	u_admin15	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin16	https://pma.studio.pathomation.com/	u_admin16	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin17	https://pma.studio.pathomation.com/	u_admin17	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin18	https://pma.studio.pathomation.com/	u_admin18	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin19	https://pma.studio.pathomation.com/	u_admin19	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm
u_admin20	https://pma.studio.pathomation.com/	u_admin20	Monday, March 27, 2022, 4:30:19 pm	Monday, March 27, 2022, 10:58:56 pm

From: <https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link: [https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=admin\\_users&rev=1648666632](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=admin_users&rev=1648666632)

Last update: **2022/03/30 21:57**

