

Server visits (administration)

This page list down all the different visits to PMA.Studio. Records can be filtered through PMA.core server url, username and date of visit. Each row consists of:

1. PMA.core server url accessed
2. PMA.Studio username used for login
3. Date of visit
4. Status of success/failure on login
5. IP address through which the user connected
6. Browser/Device used to login

The screenshot shows the Pathomation administration interface. On the left is a sidebar with navigation options: Tile Servers, Security, Panels Layout, Ribbon Layout, Global Settings, Shared Links, Conference Settings, Server Visits (highlighted), Users, Manual, Swagger Docs, and Viewer. The main area contains a search filter with fields for 'Server url' and 'Username', and a 'Search' button. Below the filter is a table with the following columns: Server URL, Username, Date, Success, IP Address, and User agent.

Server URL	Username	Date	Success	IP Address	User agent
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Wednesday, March 30th 2022, 11:23:53 pm	Success	206.84.150.33	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Wednesday, March 30th 2022, 9:53:10 pm	Success	206.84.150.33	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Wednesday, March 30th 2022, 9:29:46 pm	Success	206.84.150.33	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Wednesday, March 30th 2022, 1:14:34 am	Success	206.84.148.73	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Wednesday, March 30th 2022, 12:40:40 am	Success	206.84.148.73	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	alena	Sunday, March 27th 2022, 10:21:44 pm	Success	188.17.163.191	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Sunday, March 27th 2022, 10:21:07 pm	Success	206.84.148.73	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	b_user1	Sunday, March 27th 2022, 10:17:10 pm	Success	206.84.148.73	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36
https://devtest.pathomation.com/test/pma.core/2/	alena	Sunday, March 27th 2022, 10:15:35 pm	Success	188.17.163.191	Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/99.0.4844.84 Safari/537.36

From: <https://docs.pathomation.com/pimsls/3.1.1/> - PIMS LS 3.1.1

Permanent link: https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=admin_visits

Last update: 2023/04/06 22:59

