

Server visits (administration)

This page list down all the different visits to PMA.Studio. Records can be filtered through server url, username and date of visit. Each record consists of:

- 1- PMA.core server url accessed
- 2- PMA.Studio username used for login
- 3- Date of visit
- 4- Status of success/faliure on login
- 5- IP address through which the user connected
- 6- Browser/Device used to login

Server URL	Username	Date	Status	IP Address	User Agent
https://docs.pathomation.com/pimsls/3.1.1/	admin	Wednesday, March 30th 2022, 11:02:52 am	Success	206.84.162.30	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)
https://docs.pathomation.com/pimsls/3.1.1/	admin	Wednesday, March 30th 2022, 11:03:19 am	Success	206.84.162.30	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)
https://docs.pathomation.com/pimsls/3.1.1/	admin	Wednesday, March 30th 2022, 11:03:46 am	Success	206.84.162.30	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)
https://docs.pathomation.com/pimsls/3.1.1/	admin	Wednesday, March 30th 2022, 11:45:41 am	Success	206.84.162.30	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)
https://docs.pathomation.com/pimsls/3.1.1/	admin	Sunday, March 27th 2022, 10:24:30 am	Success	188.117.162.104	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)
https://docs.pathomation.com/pimsls/3.1.1/	admin	Sunday, March 27th 2022, 11:27:57 am	Success	206.84.162.30	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)
https://docs.pathomation.com/pimsls/3.1.1/	admin	Sunday, March 27th 2022, 11:18:30 am	Success	188.117.162.104	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/98.0.4759.25; Safari/537.36)

From:
<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:
https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=admin_visits&rev=1648665512

Last update: **2022/03/30 21:38**

