


## Version retrieval (automation)

For monitoring purposes, it is useful to determine the PMA.studio version used through a back-end API call. All of Pathomation's platform components have this option, actually.

In order to determine the installed version of PMA.studio, you can issue the call:

<https://server/pma.studio/api/json/getversioninfo>

The result is a JSON string formatted as A.B.C.D, with A, B, C being numeric. D can be either numeric or alphanumeric.



Reporting this version number to Pathomation when [asking for support through our helpdesk](#) is always helpful, too.

From:  
<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:  
[https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=auto\\_version&rev=1643105135](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=auto_version&rev=1643105135)

Last update: **2022/01/25 13:05**

