


# Other Tools

Home tab also consists of independent icons like, Synx, grid overlays, toggle colors or search panel tools.

## Synchronize Tool

The Synchronize tool contains tools for turning on synchronized viewing and some grid management tools.



## Synchronized viewing

Turning on the Sync icon when you have 2 or more cells open in your grid will trigger synchronized viewing, you can see a quick video of this below:

From:  
<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:  
[https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=other\\_icons&rev=1747060795](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=other_icons&rev=1747060795)

Last update: **2025/05/12 17:39**

