

## Whats new?

### PIMS LS 3.0.2 Patch

#### Added option to allow importing PNG and TIFF images

- Added functionality to allow PNG and TIFF images annotation import. Previously application was handling well know text(.wkt), visio pharam (.mld), indica labs (.halo) and .csv formats.

From:

<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=release\\_notes&rev=1759844001](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=release_notes&rev=1759844001)

Last update: **2025/10/07 16:33**

