

Whats new?

PIMS LS 3.0.2 Patch

Added option to allow importing PNG and TIFF images

- Added support for importing annotations from PNG and TIFF image formats.
- Previously, the application supported annotation imports only from Well-Known Text (.wkt), Visio Pharam (.mld), Indica Labs (.halo), and .csv formats.

From:

<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:

https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=release_notes&rev=1759845241

Last update: **2025/10/07 16:54**

