

## Whats new?

### PIMS LS 3.1.0

#### Features

- Added auto navigate to particular folder using link sharing. Previously slide sharing was possible through the url by manually adding parameters to it. This version provides option to generate share link for both ready only and full access control to user slides or folders.

#### Bug Fixes

- PMA.Search is no long applicable so removed it from admin new server addition/edit page.
- Fixed UX toggle issues while navigating between different slides. Now the application will retain state of toggles while navigating between slides.
- Fixed viewer version issue under About us section.

From:

<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=release\\_notes&rev=1763467155](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=release_notes&rev=1763467155)

Last update: **2025/11/18 14:59**

