


# Ribbon

We strive to make the user interface (UI) of PMA.studio (and PMA.vue) as recognizable and easy to use as possible.

The top-part of the window is a toolbar implemented as a ribbon. The ribbon is a standard application implement now that is found in both general purpose applications like Microsoft Office, and even in specialized ones like geographic information systems (GIS).



A special mentioning is warranted for the ribbon's context-menu in the left-hand corner of the ribbon: this brings up a separate menu with application-wide settings.

From:  
<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:  
<https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon&rev=1672140868>

Last update: **2022/12/27 14:34**

