


# Ribbon

We strive to make the user interface (UI) of PMA.studio (and PMA.vue) as recognizable and easy to use as possible.


The top-part of the window is a toolbar implemented as a ribbon. The ribbon is a standard application implement now that is found in both general purpose applications like Microsoft Office, and even in specialized ones like geographic information systems (GIS).



The ribbon is like a header at the top of the PMA.studio (and PMA.vue) web application (which in turn runs in a web-browser). It has tabs where all the options are categorized. Each tab has a group of options. You can use each option direct from the ribbon using the buttons.

## Context menu

A special mentioning is warranted for the [ribbon's context-menu](#) in the top left-hand corner of the ribbon: this brings up a separate menu with application-wide settings.



## Tabs

The multiple tabs, icons, groups, and the various options with their functionalities are described elsewhere in this wiki and can be consulted in two ways:

- [Tab inventory](#)
- [Command inventory](#) (useful in case you wish you [\(re\)configure your ribbon](#))

The Home tab is the default tab in Excel. It has the most frequently used options which have different groups like a clipboard, font, alignment, number, styles, cells, and editing.

From:

<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:

<https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon&rev=1746974188>

Last update: **2025/05/11 17:36**

