

# Ribbon > Home tab > Layout group > Default layout


## Button identification



## Associated command

## Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on a left and a [viewport panel](#) on the right.



## Related functionality

[Ribbon > Home tab](#)

Last update: 2023/02/08 15:44 ribbon\_home\_layout\_default [https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon\\_home\\_layout\\_default&rev=1675860240](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon_home_layout_default&rev=1675860240)

---

From: <https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link: [https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon\\_home\\_layout\\_default&rev=1675860240](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon_home_layout_default&rev=1675860240)

Last update: **2023/02/08 15:44**

