

# Ribbon > Home tab > Viewport group > Toggle slide label

## Button identification




## Associated command

viewport-toggle-label


## Function description

This function shows or hides the full path reference to the current slide in the top-right area of the ribbon

When enabled, the ribbon looks like this:



When disabled, the widget is hidden from the ribbon:



Unlike other toggle functionality, this one does not impact the viewport area

## Related functionality

[Ribbon > Home tab](#)

From: <https://docs.pathomation.com/pimsls/3.1.1/> - PIMS LS 3.1.1

Permanent link: [https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon\\_home\\_viewport\\_toggle\\_label](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon_home_viewport_toggle_label)

Last update: 2023/02/10 20:39

