

Ribbon > Home tab > Viewport group > Toggle slide label

Button identification



Associated command

viewport-toggle-label

Function description

This function shows or hides the full path reference to the current slide in the top-right area of the ribbon

When enabled, the ribbon looks like this:



When disabled, the widget is hidden from the ribbon:



Unlike other toggle functionality, this one does not impact the viewport area

Related functionality

[Ribbon > Home tab](#)

From:

<https://docs.pathomation.com/pimsls/3.1.1/> - PIMS LS 3.1.1

Permanent link:

https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon_home_viewport_toggle_label&rev=1676050515

Last update: 2023/02/10 20:35

