

# Ribbon > Home tab > Viewport group > Reset panel layout

## Button identification




## Associated command

reset - layout


## Function description

This resets the panel layout to the default settings, as defined in the [admin interface](#).

When your layout is like this e.g.:



After clicking the button the organization of the panels resets to:



The currently loaded slides are kept.

Also note that [the configuration of the panels as stipulated through the admin interface](#) can be different from the one defined by the [Default layout button](#).

## Related functionality

[Ribbon > Home tab](#)

[Default layout button](#)

[Administrative interface to define default panel configuration](#)

From:  
<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:  
[https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon\\_home\\_viewport\\_toggle\\_reset\\_panel\\_layout](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=ribbon_home_viewport_toggle_reset_panel_layout)

Last update: **2023/02/10 18:40**

