


The Ribbon

What is the ribbon?

The ribbon is the part you see at the top of PIMSLs, it houses the tools you need to use.




How is it organized?

The ribbon is broken up into two rows, similar to MS Office products.

In the top row (red box) you can see the Tabs - each one houses different functions which are explained in their own pages in the wiki.

In the picture below, you can see the Home tab is selected, indicated by the outline of the word 'Home'.

In the bottom row (blue box) you can see the tools for the Tab you're currently in.



From:

<https://docs.pathomation.com/pimsls/3.1.1/> - **PIMS LS 3.1.1**

Permanent link:

https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=the_ribbon

Last update: **2025/05/10 18:31**

