


The Viewport

This article will explain in more detail how to use the Viewport

Opening your first image

PIMS LS will always start by default with the Slide Browser panel open on the left hand side of your screen (red box).

All the files your admin has made available to you should appear in the folder tree-view. Navigate to an image and click it or drag it into the viewport (the large white space in the middle of the screen) to open it:



From: <https://docs.pathomation.com/pimsls/3.1.1/> - PIMS LS 3.1.1

Permanent link: https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=the_viewport&rev=1730991412

Last update: 2024/11/07 17:56

