

Logs

The log file section provides access to the raw log files produced by the server. This page always displays logs files available in a reverse chronological order so the latest log file is always on top. You can change the number of log files presented using the *Show Top* textbox.

There are two types of log files available:

- PMA.core log: Those are PMA.core's internal logging system with information about *operations on the server, exceptions, errors, access log* etc.
- IIS log: Those are raw IIS log files, containing all the requests made to the server along with all IIS logged info like *time taken, host, http status code* etc

Available log files		
Home > Log files		
Top 20 most recent files Show top 20		
Date	PMA Core log	IIS log
Tuesday, March 29, 2022	tile_server_log-2022-03-29.log	u_ex220329.log
	tile_server_log-2022-03-28.log	u_ex220328.log
Monday, March 28, 2022	tile_server_log-2022-03-27.log	u_ex220327.log
Sunday, March 27, 2022	tile_server_log-2022-03-26.log	u_ex220326.log
Saturday, March 26, 2022	tile_server_log-2022-03-25.log	u_ex220325.log
Friday, March 25, 2022	tile_server_log-2022-03-24.log	u_ex220324.log
Wednesday, March 23, 2022	tile_server_log-2022-03-23.log	u_ex220323.log

More insight into errors is provided via [Elmah](#).

From:
<https://docs.pathomation.com/pma.core/2.0.0/> - **PMA.core 2.x**

Permanent link:
<https://docs.pathomation.com/pma.core/2.0.0/doku.php?id=logs&rev=1648804926>

Last update: **2022/04/01 12:22**



