

User states

A PMA.core user can be in one or more states defined here

Administrator

The user is an Administrator, which has access to the Administrator UI of PMA.core as well access to all admin API features

Suspended

The user has suspended all operations, he cannot login or perform any API operation. This is as good as if a user deleted.

Locked out

The user has been locked out of the system because of too many incorrect attempts at login or by an administrator

From:

<https://docs.pathomation.com/pma.core/2.0.0/> - **PMA.core 2.x**

Permanent link:

https://docs.pathomation.com/pma.core/2.0.0/doku.php?id=user_states&rev=1648644050

Last update: **2022/03/30 15:40**

