


## Active sessions

In the session management section the currently active users are displayed. The system gives the ability to terminate a user's session.



| Home > Session Management          |                          |                    |                |                    |                             |
|------------------------------------|--------------------------|--------------------|----------------|--------------------|-----------------------------|
| Session timeout value = 20 minutes |                          |                    |                |                    |                             |
| User Name                          | Session ID               | Last Updated       | Service access | Expiration         |                             |
| admin                              | 8y1D6pO5OGQTVEJ6Hs4Z8g== | 30/1/2014 12:00:57 | True           | 30/1/2014 12:20:57 | <a href="#">✖ Terminate</a> |

The lifespan of a session is governed by the ASP.Net session time out value. This value is displayed in the installation check page. When a session is idle for more time than what is defined as the session time out value, it becomes automatically invalid and the corresponding user will need to log in again.

You may delete an active session (e.g. when it's been idle for a large period of time; perhaps the end-user forgot to log out properly), but you cannot delete your own session.

From:

<https://docs.pathomation.com/pma.core/2.0.1/> - **PMA.core 2.x**

Permanent link:

[https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=active\\_sessions&rev=1591181307](https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=active_sessions&rev=1591181307)

Last update: **2020/06/03 13:48**

