

## Form management

In the forms management section, data input forms can be designed that are subsequently assigned to virtual slides. These forms can then be filled in by the system's users while viewing slides. The data that users submit for each form are also available within this section.

From:

<https://docs.pathomation.com/pma.core/2.0.1/> - **PMA.core 2.x**

Permanent link:

[https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=form\\_management&rev=1591179198](https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=form_management&rev=1591179198)

Last update: **2020/06/03 13:13**

