



License expiration

When your license has expired or has become invalid for any other reason, you can still keep using the software. You receive a notification of this after logging in:



In addition to this administrative message, the tiles served through the tile server show as “faded”. Each tile is also provided with a notification warning, making it rather obvious and that you should [contact Pathomation](#) to resolve whatever issue is the source of the problem.




First time users

When starting the application for the first time a license will be required.


Close to expiration

When license expiration is less than a month away, a warning message appears:



Expired licenses

When the license is expired, you are notified of this right after logging in:



From:

<https://docs.pathomation.com/pma.core/2.0.1/> - **PMA.core 2.x**

Permanent link:

https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=license_expiration

Last update: **2022/03/07 14:48**

