

Server settings

In this section, various aspects of the server behavior can be configured and controlled.

Email settings

Server Settings

Home > Server settings

Email Settings

User settings

Password generation

LDAP Settings

OAuth Settings

Installation History

Use Two Factor Authentication

☐

The name of the two factor authentication

PMA.Core

Email Address

info@pathomation.com

This address will show up as the sender of the emails

Send emails using

SendGrid

SendGrid Api Key

SG.D6giWR_7HRqOXcn7pFRDaaA.JGut

Save

Clear settings

Test email

Mail log

Mail log

Mail Log

Time	Username	Description
5/25/2021, 11:59:09 AM	korinna	Mail sent with subject: PMA.control - Registration Successful
6/7/2021, 2:01:03 PM	charlotte	Mail sent with subject: Password Reminder
7/12/2021, 12:36:02 PM	charlotte	Mail sent with subject: Password Reminder
8/9/2021, 9:58:10 AM	wordpress	Mail sent with subject: Password Reminder
2/16/2022, 3:14:50 PM	veerle	Mail sent with subject: Password Reminder

Showing 1 to 5 of 5 entries

Previous

1

Next

User settings

User Settings

Home > User settings

Email Settings	User settings	Password generation	LDAP Settings	OAuth Settings	Installation History
----------------	---------------	---------------------	---------------	----------------	----------------------

Max invalid login attempts before lock out

Automatically lift lock out after (seconds)

Allow duplicate emails ☒

Force change password after days
Use 0 to disable

Password must be different from the last used
Use 0 to disable

Password generation

Password generator Settings

Home > Password settings

Email Settings	User settings	Password generation	LDAP Settings	OAuth Settings	Installation History
----------------	---------------	---------------------	---------------	----------------	----------------------

Case ☐ Lower ☐ Upper ☒ Mixed ☐ Any

Length

Include letters ☒

Include numbers ☒

Include special characters ☐

LDAP settings

LDAP Settings

Home > LDAP settings

Email Settings	User settings	Password generation	LDAP Settings	OAuth Settings	Installation History
----------------	---------------	---------------------	---------------	----------------	----------------------

Use LDAP ☐

LDAP Server	LDAP Root	LDAP Username	LDAP Password
-------------	-----------	---------------	---------------

Add new LDAP server

Add new LDAP Server

LDAP Server

LDAP distinguished name

Username

Password

Save

[Back to List](#)

OAuth settings

OAuth Settings

Home > OAuth settings

Email Settings | User settings | Password generation | LDAP Settings | **OAuth Settings** | Installation History

Callback URL: <https://localhost:3000/auth/callback?next=/#/Login/CustomCallback>

OAuth Name	OAuth Grant Type	OAuth Public Key	Auth url	Token url	Client	
Identity Server	sub		https://localhost:3000/auth/callback?next=/#/Login/CustomCallback	https://localhost:3000/auth/callback?next=/#/Login/CustomCallback	Public	 

 Add new OAuth server

Add new OAuth server

Add new OAuth Server

Name		
Public Key		The Public Key field is required
User Id claim		The User Id claim field is required
Auth Uri		
Token Uri		
Client Id		
Client Secret		Leave empty if no client secret is required for code to token exchange
Callback Uri	https://localhost:44300/authorize-redirects/1234567890/OAuthCallback	

[Save](#) [Back to List](#)

Installation history

Installation History									
Home > Installation History									
Email Settings	User settings	Planned generation	LDAP Settings	OAuth Settings	Installation History				
Date	Version	MY Version	PHP 2.0 Version	License	Downgrade	Repaired code	License change date		
2025/05/02 5:17:38 PM	2.0.2.20251	1.3.8	2.16.0	☐	dan	code repaired by the system			
2025/05/02 9:23:54 AM	2.0.1.20250	1.3.8	2.16.0	☐		code repaired by the system			
2025/05/02 9:29:17 AM	2.0.1.1602	1.3.8	2.16.0	☐		code repaired by the system			

From:

<https://docs.pathomation.com/pma.core/2.0.1/> - **PMA.core 2.x**

Permanent link:

https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=settings_server&rev=1648486241

Last update: **2022/03/28 19:50**

