





Submitted data

The data submitted by each user for each form can be obtained in the “Submitted data” tab in the formats:

- Unordered List ItemHTML
- Unordered List ItemCSV
- Unordered List ItemXML
- Unordered List ItemARFF (Atribute-Relational file format)



The screenshot shows a web interface titled "Forms". Below the title is a breadcrumb trail: "Home > Forms > Submitted data". There are two tabs: "Submitted data" (active) and another tab with a grid icon. Below the tabs is a table for the "NSCLC form". The table has columns for "Version", "Results", and four columns with icons representing different data formats: a calendar icon, a document icon, a CSV icon, and an "Arff" link.

NSCLC form					
Version	Results				
1	2170				Arff

From:

<https://docs.pathomation.com/pma.core/2.0.1/> - PMA.core 2.x

Permanent link:

https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=submitted_data&rev=1591180417

Last update: 2020/06/03 13:33

