

Submitted data

The data submitted by each user for each form can be obtained in the “Submitted data” tab in the formats:

- HTML
- CSV
- XML
- ARFF (Attribute-Relational file format)

Forms						
Home > Forms > Submitted data						
Submitted data						
NSCLC form						
Version		Results				
1		2170				 Arff

From:
<https://docs.pathomation.com/pma.core/2.0.1/> - **PMA.core 2.x**

Permanent link:
https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=submitted_data&rev=1646127980

Last update: **2022/03/01 12:46**

