

Audit trail

The audit trail section displays a log of modification events that involve the current user.

1230 x 471 px (displayed at 854 x 327)

Audit trail for user: test

Home > Users > Audit trail

Summary Edit Activity log **Audit trail** ACL

EntityId	EventId	ActionType	EventDate	UserLogin	UserHost	Administrator	CanAnnotate	Description	Email	FirstName	LastName
0	174	Update	4/15/2019, 4:34:15 PM	admin	85.73.248.94	False	True	Test user	test123@pathomation.com	Test	test
0	173	Create	4/15/2019, 4:34:00 PM	admin	85.73.248.94	False	False	Test user	test123@pathomation.com	Test	test

Showing 1 to 2 of 2 entries

Previous 1 Next

From:
<https://docs.pathomation.com/pma.core/2.0.2/> - **PMA.core 2.x**

Permanent link:
https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=audit_trail&rev=1591175016

Last update: **2020/06/03 12:03**

